




>>Entry Name: 4,4-Diphenylpiperidine



CAS Registry Number: 34273-01-3

Molecular Formula: C<sub>17</sub>H<sub>19</sub>N

Molecular Weight: 237.344

Accurate Mass: 237.151749

Percentage Composition: C 86.03%; H 8.07%; N 5.90%

Melting Point: Mp 72-74°

Boiling Point: Bp<sub>0.005</sub> 129-131°

>>Derivative: Hydrochloride

Melting Point: Mp 300-305°

>>Derivative: *N*-tert-Butyl

Synonym(s): 1-*tert*-Butyl-4,4-diphenylpiperidine. **Budipine**, INN. BY 701

CAS Registry Number: 57982-78-2

Type of Compound Code(s): XA2050 XA1450

Molecular Formula: C<sub>21</sub>H<sub>27</sub>N

Molecular Weight: 293.451

Accurate Mass: 293.214349

Percentage Composition: C 85.95%; H 9.27%; N 4.77%

Use/Importance: NMDA antagonist. Antiparkinsonian agent. Antidepressant

Development Status: In Phase III clinical trials (1993)

Melting Point: Mp 108-109°

Boiling Point: Bp<sub>0.05</sub> 160-165°

Calculated Log P Data: Log P 5.6 (calc)

**Hazard & Toxicity:** LD<sub>50</sub> (rat, orl) 165 mg/kg. LD<sub>50</sub> (rat, ivn) 28 mg/kg

>>Derivative: *N*-tert-Butyl; hydrochloride

Synonym(s): Parkinsan

CAS Registry Number: 63661-61-0

Melting Point: Mp 306-309 dec.°

#### References:

- Menge, H.G. *et al.*, *Arzneim.-Forsch.*, 1982, **32**, 85, (*pharmacol*)  
Schaefer, H. *et al.*, *Arzneim.-Forsch.*, 1984, **34**, 233, (*synth, pharmacol*)  
Gerstenbrand, F. *et al.*, *Clinical Exper. Budipine Parkinson Ther. [Proc. Int. Workshop Parkinson's Dis.]*, Springer-Verlag, Berlin, 1985, (*book*)  
Caputo, O. *et al.*, *Eur. J. Drug Metab. Pharmacokinet.*, 1991, **16**, 113, (*metab*)  
Klockgether, T. *et al.*, *J. Neural Transm.*, 1993, **5**, 101, (*pharmacol*)  
Martindale, *The Extra Pharmacopoeia*, 30th edn., Pharmaceutical Press, 1993, 840  
Jackisch, R. *et al.*, *Eur. J. Pharmacol.*, 1994, **264**, 207, (*pharmacol*)  
Klumpp, D.A. *et al.*, *J.O.C.*, 1999, **64**, 6702-6705, (*synth, pmr, cmr*)  
Lewis, R.J., *Sax's Dangerous Properties of Industrial Materials*, 8th edn., Van Nostrand Reinhold, 1992, BOM530