



>>Entry Name: Tetraiodomercurate(2-)

CAS Registry Number: 14876-64-3

[HgI₄]²⁽⁻⁾

Molecular Formula: HgI₄²⁽⁻⁾

Molecular Weight: 708.208

Accurate Mass: 709.588509

Percentage Composition: Hg 28.32%; I 71.68%

>>Derivative: Di-K salt

Synonym(s): Dipotassium tetraiodomercurate. Potassium mercuric iodide

CAS Registry Number: 7783-33-7

Molecular Formula: HgI₄K₂

Molecular Weight: 786.405

Accurate Mass: 787.515923

Percentage Composition: Hg 25.51%; I 64.55%; K 9.94%

Use/Importance: Dil. soln. in NaOH is *Nessler's reagent*, used for NH₄⁽⁺⁾ analysis

Biological Use/Importance: Antiseptic

Aldrich: 34514-8

Fluka: 72190

>>Derivative: Di-K salt, bis(dibenzo-18-crown-6) complex

Physical Description: Cryst.

>>Derivative: Bis-Me₄N salt

CAS Registry Number: 33615-18-8

Physical Description: Cryst.

References:

Mok, K.F., *Inorg. Chim. Acta*, 1988, **142**, 27, (synth, ir, pmr)

Pakhomov, V.I. *et al.*, *Zh. Neorg. Khim.*, 1991, **36**, 150, (synth, cryst struct)

Lewis, R.J., *Sax's Dangerous Properties of Industrial Materials*, 8th edn., Van Nostrand Reinhold, 1992, NCP500