




>>Entry Name: Xylenol orange

Synonym(s): *N,N'*-[3*H*-2,1-Benzoxathiol-3-ylidenebis[(6-hydroxy-5-methyl-3,1-phenylene)methylene]]bis[*N*-(carboxymethyl)glycine] *S,S*-dioxide, 9Cl. 3,3'-Bis[*N,N*-di(carboxymethyl)aminomethyl]-*o*-cresolsulfonephthalein. *o*-Cresolphthalexone *S*



CAS Registry Number: 1611-35-4

Related CAS Registry Numbers(s): 3618-43-7

Molecular Formula: C₃₁H₃₂N₂O₁₃S

Molecular Weight: 672.665

Accurate Mass: 672.162513

Percentage Composition: C 55.35%; H 4.79%; N 4.16%; O 30.92%; S 4.77%

Use/Importance: Used as 0.05% aq. soln. for spectrophotometric detn. of Bi ($\lambda_{\max}$ 545 nm, ϵ 24000), Zr ($\lambda_{\max}$ 535 nm, ϵ 35000), Hg, Nb, Sc, Y, Zn, Al, Th, Ga, U, V; as an indicator in EDTA titrations. Mono- and tetrasodium salts are also used

Physical Description: Dark red cryst. powder

Melting Point: Mp 286 dec.°

Solubility: Sol. H₂O

pKa Value: p*K*_{a5} 2.58; p*K*_{a7} 6.4; p*K*_{a8} 10.5 (25°)

>>Derivative: Mono-Na salt

CAS Registry Number: 63721-83-5

Physical Description: Dark red cryst. powder

Melting Point: Mp 210 dec.°

Other Data: $\lambda_{\max}$ 580 nm

Fluka: 95615

References:

Körbl, J. *et al.*, *Chem. Anal. (Warsaw)*, 1956, **45**, 102, (*indicator*)

Cheng, K.L., *Talanta*, 1959, **2**, 61; 266; **3**, 81, (*detn, Zr, Hf*)

Budešinsky, B., *Coll. Czech. Chem. Comm.*, 1962, **27**, 226, (*detn, Th*)

Cabrera-Martin, A. *et al.*, *Talanta*, 1969, **16**, 1023, (*detn, Hg*)

Leong, C.L., *Anal. Chem.*, 1973, **45**, 21, (*detn, U*)

Snell, F.D., *Photometric and Fluorimetric Methods of Analysis, Metals, Part I*, John Wiley, New York, 1978, 523; *Metals, Part II*, 1062; 1167; 1923

Onishi, H., *Photometric Determination of Traces of Metals, Part IIa, Individual Metals. Aluminium to Lithium*, John Wiley, New York, 4th Ed., 1986, 787; *Part IIb, Individual Metals. Magnesium to Zinc*, John Wiley, New York, 4th Ed., 1989, 164; 357; 767

Marczenko, Z., *Separation and Spectrophotometric Determination of Elements*, Horwood, Chichester, 1986, 172; 501; 648, (*use*)