



>>Entry Name: 3,5-Diaminobenzoic acid, 9CI

CAS Registry Number: 535-87-5

Molecular Formula: $C_7H_8N_2O_2$

Molecular Weight: 152.152

Accurate Mass: 152.058578

Percentage Composition: C 55.26%; H 5.30%; N 18.41%; O 21.03%

Use/Importance: Used in the fluorimetric determination of DNA. Monomer in the manuf. of polyamides used as reverse osmosis membranes

Physical Description: Needles

Melting Point: Mp 228° (slow heat). Mp 240° (rapid heat)

Solubility: Sl. sol. H_2O ; sol. oxygenated solvs.

pKa Value: pK_a 5.3 (25°)

Other Data: Forms a monohydrate

Aldrich: D1280-5

Fluka: 32770

Sigma: D8251

>>Derivative: Hydrochloride (1:2)

CAS Registry Number: 618-56-4

Melting Point: Mp 300°

Aldrich: 11383-2

Fluka: 32780

Sigma: D1891

>>Derivative: Et ester

CAS Registry Number: 1949-51-5

Molecular Formula: $C_9H_{12}N_2O_2$

Molecular Weight: 180.206

Accurate Mass: 180.089878

Percentage Composition: C 59.99%; H 6.71%; N 15.55%; O 17.76%

Melting Point: Mp 84°

References:

Aldrich Library of ^{13}C and 1H FT NMR Spectra, 1992, **2**, 1138A; 1138B, (*nmr*)

Aldrich Library of FT-IR Spectra, 1st edn., 1985, **2**, 227B; 227C, (*ir*)

Aldrich Library of IR Spectra, 2nd Ed., 856D; 856E, (*ir*)

Aldrich Library of NMR Spectra, **6**, 174D; 175A, (*pmr*)

Kailan, A. *et al.*, *Monatsh. Chem.*, 1903, **56**, 407

Robertson, F.W. *et al.*, *Anal. Biochem.*, 1971, **41**, 477

U.S.S.R. Pat., 1971, 311 903; *CA*, **75**, 140498s, (*synth*)

Hinderer, H.E. *et al.*, *Appl. Polym. Symp.*, 1973, **21**, 1, (*polyamides*)

Lepri, L. *et al.*, *J. Chromatogr.*, 1974, **88**, 331, (*tlc*)

Sax, N.I., *Dangerous Properties of Industrial Materials, 5th edn.*, Van Nostrand Reinhold, 1979, 541